MASAJEADOR FACIAL

MANUAL DE INSTRUCCIONES TH – BY23

ATENCIÓN.

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INSTRUCTIONS MANUAL TH – BY23

ATTENTION.

- Please read this manual carefully before using the device for the first time.
- Please keep this instruction manual for future reference.
- If you lend or sell the device to third parties, also include the manual.

Dear Customer:

We appreciate the trust placed in us when purchasing this THULOS® product. Our main goal is to offer excellent items, maintaining the strictest quality controls and complying with all current safety regulations, with the aim of presenting a perfect and safe product for home use.

GENERAL SAFETY INSTRUCTIONS.

- Failure to comply with these safety regulations will be grounds for waiving the warranty.
- Never pull on the cable to unplug the charging connector, hold and pull on the connector.
- Do not place the appliance near hot elements such as stoves, ovens, etc.
- Never leave the device unattended while it is plugged in and charging.
- This appliance is intended for home and occasional use only.
- Do not use outdoors, this device is designed for indoor use.
- The device can be used by children from 8 years of age and by people with reduced physical, sensory or mental capabilities, or who do not have experience or knowledge, if they are supervised or instructed in the use of the appliance safely and they understand the risks involved.
- Children cannot play with the device or clean it.
- The father, mother or guardian will be fully responsible for the child using the device.
- Never immerse the device, charging cable or adapter in water or any other liquid.
- Unplug the device immediately after fully charging.
- Do not use accessories not recommended by the manufacturer, they may cause damage.
- If you observe any anomaly in the device, contact the official technical service.
- Never try to repair the device yourself, always contact your technical service.
- Do not place the device on surfaces that may be hot.
- Do not touch the device with wet hands or body or barefoot while it is charging.
- Always keep the device clean after each use.

PRECAUTIONS.

- Do not use this device if you have had any type of plastic surgery, if you have subcutaneous metal implants or if you are pregnant.
- Only to use on the face and neck, avoiding acne and redness.
- People with allergic or atopic skin should not use this device daily, but can be used once a week.
- Each mode can be used only once a day, but we recommend that you successively use one mode each day.

****BATTERY****

- If the battery runs out during use, charge immediately.
- Always keep the battery charged even if you do not use the device.
- If not using, charge the battery once a month.

PARTS OF THE APPLIANCE.

1. Head.
2. Indicator light.
3. On/Off and mode selection.
4. Light (red, blue or green).
5. Function indicator.
6. Charging connector.

OPERATION.

1. Charge the device fully (charging time will depend on the adapter used).
2. Clean your face before applying any treatment.
3. To turn on the device, press and hold the button until it turns on.
(Won't turn on while plugged in charging).
4. Select the desired light and massage mode by short pressing the button.
5. Place the massage head on the area to be treated and the massage will begin.
6. To turn off the device, press and hold the button until it turns off.

CLEANING AND MAINTENANCE.

- Never clean the device while it is charging.
- Never submerge the device in water. Clean with a slightly damp cloth.
- You can moisten the cloth in a soapy solution of water and neutral soap.
- Dry with a soft, clean cloth or absorbent kitchen paper.
- Never use scouring pads or abrasive cleaning products.
- Do not put the appliance in the dishwasher.
- Store in a safe and dry place, free of dust. If possible in its original box.

CLEANER MODE (CLEAN Option, green light):

This mode can effectively clean the pores. Combined with green light it can improve cellular oxygenation and facilitate microcirculation.

Turn on the device and press the button until the CLEAN option is activated with a green light. Apply your cleanser or makeup remover lotion to your skin (do not apply directly to the apparatus). Place the head on your skin and gently move it from the bottom up.

TREATMENT MODE (EMS Option, red light):

This mode can improve the appearance of the skin and helps facial care products absorb better. Combined with red light it helps revitalize the skin and reduce the appearance of wrinkles.

Turn on the device and press the button until the EMS option is activated with a red light.

Apply your facial toner or lotion (do not apply directly to the device).

Place the head on your skin and gently move it from the bottom up.

HEAT MODE (HOT option, blue light):

This mode helps relax the skin. Combined with blue light it helps to rejuvenate and to redefine the skin.


Turn on the device and press the button until the HOT option is activated with blue light.

Apply your mask to your face and neck (do not apply directly to the device).

Place the head on your skin and gently move it from the bottom up.

TIPS.

We recommend massaging following these schemes:

SPECIFICATIONS.

- Compatible adapter: 5V up to 2 Amps maximum. (Not included).
- Usb-A to MicroUsb charging cable (Included).
- Charging time: will depend on the amperage of the adapter used.
- Massage vibration power: 2.5 W.
- Battery: 3.7V – 500mAh.
- 3 Facial Treatment Modes.
- 7 Levels of massage intensity.

WARRANTY CONDITIONS.

- This device is covered against any breakdown as long as it has been caused due to failure of the manufacturer's own electrical system and is within the period warranty.
- The warranty does not cover knocks, breakdowns due to misuse of the device or any other negligence by the user. The manufacturer reserves the right to void the warranty in the moment any negligence on the part of the user is observed.
- The documents that validate the guarantee can be the invoice or purchase receipt.

where the brand and model of the device or its warranty card are reflected duly completed and sealed by the seller.

- During the warranty period, the store where you purchase the device is responsible for manage and process the guarantee in case of any breakdown.

You must take the device with its documentation to the store where you purchased it.

ENVIRONMENTAL PROTECTION STANDARDS.

This appliance should not be disposed of with household waste at the end of its useful life, but should be offered to a recycling center for household electrical and electronic equipment. This symbol on the appliance, instruction manual and packaging alerts you to this important aspect.

The materials used in this device can be recycled. By recycling household appliances, you contribute to promoting environmental protection. Please consult your local authorities for information about the collection point.

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