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Dear Customer:

We thank you very much for having purchased this THULOS® product. We kindly ask you to read and save this instruction manual. Please read these operating instructions carefully and keep them in a safe place. If you pass the unit on to someone else, please give them these operating instructions as well.


WARNING

- Before use, talk to your doctor if you have any medical issues or concerns.
- This product is not for use by children.
- Do not wear necklaces, earrings or any facial jewelry while using this product.
- Necklaces, earrings and other types of jewelry should be removed before using this product.
- This product is not for use on areas of the body that have undergone any sort of surgery or have any scars.
- If you experience pain while using this product, discontinue use immediately.
- This product is not waterproof and should not be used in the shower or any situation that involves water.
- Do not store this product in the bathroom or a room that has high moisture content.
- Do not operate in the presence of explosives and/or other flammable fumes or liquids.
- Do not store or use this product near sources of heat or fire.

WARNING

KEEP BATTERIES AWAY FROM CHILDREN HARMFUL IF SWALLOWED

PRODUCT COMPONENTS

SMALL HEAD

Gentle suction, perfect for everyday use. Good for use on sensitive skin and nose area.

LARGE HEAD

Strong suction, perfect for blackheads. Note: Not for use around the eye area.

OVAL HEAD

Oval opening for getting into corners of nose and mouth, and for sensitive skin around eyes.

SONIC HEAD

Helps exfoliate dead skin while pulling dirt & oil from skin.

Note:

For Dry Skin – Once a week.

For Oily/Mixed Skin – Twice a week.

No more than 5 minutes a day.

BATTERY INSTALLATION


Figure A

3. Open the battery cover (Figure A).

Figure B

4. Insert 2*AA batteries (not included) into the battery compartment. Make sure to line up the (+) and (-) markings on the battery with the (+) and (-) markings on the inside of the battery compartment (Figure B).

HOW TO CHANGE THE SUCTION HEADS


Figure C


Figure D

3. Place one hand on the suction head and your other hand on the unit (Figure C)

4. Firmly pull off the suction head and reinstall your desired suction head in the same manner (Figure D)

HOW TO USE


Figure E


Figure F

- 3 Start by opening up your pores by washing your face with warm water, steam or a hot wet towel (Figure E). Make sure to dry face and hands completely before moving onto the next step.

- 4 Select your desired suction head (see "How to change the suction heads"). The device is turned on by sliding the ON/OFF switch up. The middle position is Low suction and the top position is High suction (Figure F).



Figure G

3 Apply the Pore Cleaner to your face (Figure G) slowly move the device around your face, not keeping it on one place for more than 2 seconds. For best results, focus on your problem areas. NOTE: DO NOT move back and forth in same direction.

HELPFUL TIPS

- After use the Pore Cleaner make sure to apply moisturizer to the treated area.
- Before and after use, do not use any facial scrubs or other exfoliating products.
- Do not apply any downward force during use, the suction capabilities of the device are enough to clean your pores efficiently.
- Slight redness is normal and will disappear after approximately 30 minutes.

HOW TO CLEAN

After each use, clean the suction head by running under warm water (Figure H) allow the suction head to air dry. Do not attempt to clean the main unit by running under water. The clean the main unit, wipe down a soft damp cloth.

WARRANTY CONDITIONS

1. This appliance is covered against any factory defects.
2. The device is changed for a new one as long as it is within the change period and is returned with all accessories and packaging in good condition.
3. The guarantee will be valid as long as the appliance is accompanied by your invoice or purchase ticket where the exact device and model, or your warranty card duly completed and sealed, is reflected.
4. THE APPLIANCE WILL BE EXEMPT FROM WARRANTY WHEN:
 - 4.1. The appliance has signs of knocks and / or falls.
 - 4.2. The appliance shows signs of misuse.
 - 4.3. The appliance has signs of overuse.
 - 4.4. The appliance is improperly handled.
 - 4.5. No supporting document. (Card, ticket or invoice).
 - 4.6. The appliance has not been used in accordance with the instructions for use.

GUIDELINES FOR PROTECTION OF THE ENVIRONMENT

This appliance should not be put in to the domestic garbage at the end of its durability, but must be offered at a central point for there cycling of electric and electronic domestic appliances. This symbol on the appliance, instruction manual and pack again put your attention to this import and tissue.

The materials used in this appliance can be recycled. By recycling of used domestic appliances, you contribute important push to the protection of our environment. Ask your local authorities for information regarding the point of recollection.

LIMPIA POROS FACIAL
Pore Cleaning Device
TH-BY01

EAN: 8436588770981

www.thulos.com

thulos@thulos.com

Servicio Técnico / Technical Service:

sat@grupothulos.com