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MANUAL DE INSTRUCCIONES **TH-BY30**

Atención.

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INSTRUCTION MANUAL **TH-BY30**

Attention.

- Read this instruction manual carefully before using the device.
- Keep the manual for future reference.
- If you give or sell the device to others, it must be accompanied by the manual.
- Failure to comply with the instructions in this manual will void the warranty.

GENERAL SAFETY INSTRUCTIONS.

- Failure to comply with these safety rules will void the warranty.
- If the charging cable is damaged, it must be replaced by the manufacturer to avoid hazards.
- Do not place the device near hot items such as stoves, ovens, etc.
- This device is designed only for domestic and occasional use.
- The device can be used by children aged 8 years and older, and by people with reduced physical, sensory, or mental capabilities, or those without experience or knowledge, if supervised or instructed on its safe use and they understand the risks involved.
- The parent, guardian, or caregiver will be fully responsible for ensuring the child uses the device properly.
- Never submerge the cable, plug, or device in water or any other liquid.
- Unplug the device immediately after charging the battery.
- Do not use accessories not recommended by the manufacturer, as they may cause fires, electric shocks, or injuries.
- If you notice any abnormalities in the device, contact official customer service.
- Do not place the device on surfaces that may be hot.

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BATTERY WARNINGS

- Charging must be supervised by an adult.
- Do not submerge the item or the batteries/battery packs in water.
- Do not expose the batteries/battery packs/accumulators to heat or fire.
- If the liquid leaking from the batteries/battery packs/accumulators comes into contact with clothing or skin, wash the area immediately with clean water and seek medical attention right away.
- Stop using the batteries/battery packs/accumulators immediately if you notice any abnormalities.
- Avoid contact between the batteries/battery packs/accumulators and metallic objects such as keys or coins.
- The batteries/battery packs/accumulators can explode or release toxic substances. Risk of fire or burns. Do not open, crush, modify, disassemble, puncture, heat, or burn them.
- A battery pack should be understood as a single unit that should not be broken, opened, or handled under any circumstances.
- Do not put batteries/battery packs/accumulators in your mouth. If swallowed, contact your doctor or the local poison control center immediately.
- Do not cause short circuits to avoid burns or fires.
- If the batteries/battery packs/accumulators are swollen, stop using them and dispose of them properly.

Danger: Failure to comply with the above warnings could result in the risk of electric shock, leakage/spillage, overheating, explosion, or fire, and may cause serious personal injury and/or property damage.

MEDICAL WARNINGS

Do not use the massage gun without approval from your doctor if any of the following conditions apply:

- Recent surgeries of any kind: skin diseases and/or injuries caused by vibrating massage guns.
- Pregnancy.
- Diabetes/neuropathy or retinal damage.
- Pacemakers.
- Wounds or injuries.

ADVERSE REACTIONS

Skin irritation and mild bruising may occur if the massage gun is used in one localized area for more than 60 seconds or on bony areas. Adverse reactions may also occur in people recovering from surgery, with skin conditions, neurological disorders, wounds, bone fractures, and/or contusions. As a precaution to avoid such reactions, it is recommended that you consult your doctor before using the massage gun.

PARTS OF THE MASSAGE GUN

1. Battery level / massage head frequency LED, massage frequency adjustment.
2. Multifunction button: power on/off.
3. USB-C charging connector.
4. Spherical head.
5. U-shaped head.
6. Flat head.
7. Small spherical head.

TYPES OF HEADS

Spherical Head (4)

This head is designed to be used on the hands, lower back, glutes, thighs, and calf muscles. It helps accelerate muscle recovery, relaxes the muscles, and promotes the release of lactic acid.

U-Shaped Head (5)

The U-shaped head is specifically designed for the cervical area and the Achilles tendon. In just 30-60 seconds, it stimulates muscle groups, helping to relieve tension.

Flat Head (6)

This head provides a gentle and even percussion massage on various muscle areas of the body. It offers effective muscle relaxation.

Small Spherical Head (7)

Ideal for massaging the spinal area and large muscle groups, such as the thigh muscles.

HOW TO USE

1. Press and hold the multifunction button (2) for 2 seconds to turn the device on or off. When the device is on, the first LED on the left (1) will light up, indicating the minimum vibration frequency of the heads.
2. With the device on, press the multifunction button (2) to change the pulse frequency of the head. Each time the button is pressed, the number of LEDs lit (1) will increase until it reaches the maximum frequency (the device has 6 different frequencies). The next time you press it, it will return to the minimum frequency.
3. With the device off, pressing the multifunction button (1) will check the battery level, which will be 100% when all the LEDs (1) are lit.
4. To replace the brush heads, simply pull them off firmly. To reposition them, align the brush head with the vibrating part of the gun, and apply gentle pressure until the head's body makes contact with the vibrating element.

BATTERY RECHARGE

1. A 5V-1A USB adapter (not included) is required for charging.
2. Before the first use, perform a 3-hour charging cycle. Connect the charging cable to the 5V adapter to begin charging.
3. LED Indicator (1): During charging, all 6 lights will flash. When the battery is fully charged, all 6 lights will remain on. During use, the lights will show the remaining battery level. Recharge the device when only one LED is lit.

WARRANTY CONDITIONS.

- This device is covered for any electrical malfunction, provided the malfunction was caused by a failure in the manufacturer's own electrical system.
- The warranty does not cover damage from impacts, malfunctions caused by misuse of the device, or any malfunction resulting from user negligence.
- During the warranty period, in case of any malfunction, you must take the device to the store where it was purchased, along with your purchase receipt or invoice.

ENVIRONMENTAL PROTECTION RULES.

This device should not be disposed of with household waste when it reaches the end of its lifespan. Instead, it should be taken to a recycling center for electrical and electronic devices. This symbol on the device, instruction manual, and packaging alerts you to this important matter.

The materials used in this device can be recycled. By recycling electronic appliances, you contribute to promoting environmental protection. Please consult your local authorities for information on the nearest collection point.

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