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TH-BD121

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TH-BD121

BATTERY.

- Before the first use, put the batteries.
- Always use alkaline high quality batteries.
- If the scale is not used for long periods, it is advisable to remove the batteries.
- Keep batteries out of the reach of children.
- Store batteries in areas away from heat sources.
- Never mix different types of batteries such as new and old.
- If the scale does not work properly, replace the batteries with new ones.

MONITORING THROUGH APP

- Install the OKOK-International App, available for Ios and Android.
- Keep the scale on.
- Turn on the bluetooth of your mobile device.
- Enter the App and in settings search for devices and then search for the scale.
- Once the scale is linked, you can now keep track of your weight.
- Every time the scale registers a weight, it will send it to the app.
- Put all the data that the app asks for and you can keep

track of your weight.

LOCATION.

For a correct operation of the balance, it must be placed on a table firm, smooth and level. Otherwise, the weight can lead to errors.

OPERATION.

- The scale turns on automatically when it is climbed on top.
- First put one foot on it until the scale turns on.
- Wait for "0" to appear on the display.
- Once "0" is zero, you can climb on top to weigh yourself.
- You must put both feet centered, within the surface of the scale.
Any other position of the feet can give wrong weight.
- The scale will calculate your weight and the display will flash.
- Once the display is fixed, that will be its weight.
- The scale turns itself off automatically.

SPECIFICATIONS.

- Size 300x300mm
- KG / LB weighing unit
- LCD Screen.
- Material Tempered Glass 6mm thick + ABS type plastic
- Batteries: AAA Type (x2) included.
- Auto off.
- Auto Reset to Zero "0".
- Capacity / Maximum Weight 180kg
- Minimum Weight 100gr

www.thulos.com
Information: thulos@thulos.com
Technical Service: sat@grupothulos.com